## Important Notes

The load-bearing capacity of the intramedullary lengthening nail FITBONE® (TAA, TAM, SAA and FSA/TSA) is – especially during the extension phase and early healing stage – limited. Please do not put more weight-stress on your leg than recommended by your doctor. Unwanted overloads e.g. falls or stumbles should be avoided in each case. Every case of unwanted overload nevertheless has to be reported to the doctor immediately.

Any interruption of distraction may cause premature bone consolidation. The FITBONE® is not a permanent implant. It must be removed after completed distraction and healing.

The owner of this Implant ID carries an artificial implant made of metal, electronic components and plastic material. Metal detectors may respond to it. This document does not release the implant carrier from being body searched in case of a positive detector signal. Strong image artefacts can be expected around the implant area.

If injections need to be made around the implant area, please inform the doctor that this may damage the implant (e.g. the receiver).

Unsafe in Magnetic Resonance Imaging Environments. The patient implanted with the FITBONE® should not enter the MRI unit. Any form of an electrical therapy in which current is conducted through the body of the

patient must be avoided at the concerned limb; likewise therapeutic ultrasonics.

The manufacturer does not know of any reports on reactions of the electromechanical intramedullary nail due to high-energy ionizing radiation, with the exception of diagnostic X-rays. Such treatment should be omitted in any case during the distraction phase.

Using defibrillators can have an adverse effect on the intramedullary lengthening nail. That is why such applications should be avoided if possible.

If it is compulsory to apply one of the above mentioned diagnostic and therapeutic procedures, we strongly recommend to contact the doctor or the manufacturer of the implant.

All cases mentioned require careful monitoring of the implant functions.

**For further safety information please consult the Instruction for Use FITBONE® Control Set.**

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