BonAlive[®]

granules

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
Literatur

- In vitro* antibiofilm activity of bioactive glass S53P4. Drago L, Vassena C, Fenu S, De Vecchi E, Signori V, De Francesco R, Romanò CL. *Future Microbiol.* 2014;9(5):593–601.
- Bioactive glass BAG-S53P4 for the adjunctive treatment of chronic osteomyelitis of the long bones: an *in vitro* and prospective clinical study. Drago L, Romanò D, De Vecchi E, Vassena C, Logoluso N, Mattina R, Romanò CL. *BMC Infectious Diseases* 2013;13:584.
- Antibacterial effects and dissolution behavior of six bioactive glasses. Zhang D, Leppäranta O, Munukka E, Ylänen H, Viljanen MK, Eerola E, et al. *J Biomed Mater Res.* 2010;93A(2):475–483.
- Bactericidal effects of bioactive glasses on clinically important aerobic bacteria. Munukka E, Leppäranta O, Korkeamäki M, Vaahio M, Peltola T, Zhang D, et al. *J Mater Sci: Mater Med.* 2008;19(1):27–32.
- Antibacterial effect of bioactive glasses on clinically important anaerobic bacteria *in vitro*. Leppäranta O, Vaahio M, Peltola T, Zhang D, Hupa L, Ylänen H, et al. *J Mater Sci: Mater Med.* 2008;19(2):547–551.
- In situ* pH within particle beds of bioactive glasses. Zhang D, Hupa M, Hupa L. *Acta Biomaterialia.* 2008;4(5):1498–1505.
- Factors controlling antibacterial properties of bioactive glasses. Zhang D, Munukka E, Hupa L, Ylänen H, Viljanen MK, Hupa M. *Key Engineering Materials.* 2007;330–332:173–176.
- Comparison of antibacterial effect on three bioactive glasses. Zhang D, Munukka E, Leppäranta O, Hupa L, Ylänen H, Salonen J, et al. *Key Engineering Materials.* 2006;309–311:345–348.
- Interactions between the bioactive glass S53P4 and the atrophic rhinitis-associated microorganism *Klebsiella ozaenae*. Stoor P, Söderling E, Grenman R. *J Biomed Mater Res.* 1999;48(6):869–874.
- Antibacterial effects of a bioactive glass paste on oral micro-organisms. Stoor P, Söderling E, Salonen JI. *Acta Odontol Scand.* 1998;56(3):161–165.
- Interactions between the frontal sinusitis-associated pathogen *Haemophilus Influenzae* and the bioactive glass S53P4. Stoor P, Söderling E, Andersson OH, Yli-Urpo A. *Bioceramics.* 1995;8:253–258.

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Literatur

- Genetic design of bioactive glass. Hench L. J Eur Cer Soc. 2009;29:1257-1265.
Molecular basis for action of bioactive glasses as bone graft substitute. Välimäki VV, Aro HT. Scand J Surg. 2006;95(2):95-102.
Intact surface of bioactive glass S53P4 is resistant to osteoclastic activity. Wilson T, Parikka V, Holmbom J, Ylänen H, Penttinen R. J Biomed Mater Res. 2005;77A(1):67-74.
Granule size and composition of bioactive glasses affect osteoconduction in rabbit. Lindfors NC, Aho AJ. J Mater Sci: Mater Med. 2003;14(4):265-372.
Histomorphometric and molecular biologic comparison of bioactive glass granules and autogenous bone grafts in augmentation of bone defect healing. Virolainen P, Heikkilä J, Yli-Urpo A, Vuorio E, Aro HT. J Biomed Mater Res. 1997;35(1):9-17.
Long term behaviour of bioactive glass cone and granules in rabbit bone. Heikkilä JT, Salonen H, Yli-Urpo A, Aho AJ. Bioceramics. 1996;9:123-126.
Bone formation in rabbit cancellous bone defects filled with bioactive glass granules. Heikkilä JT, Aho HJ, Yli-Urpo A, Happonen R, Aho AJ. Acta Orthopaedica. 1995;66(5):463-467.

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

Literatur


Molecular basis for action of bioactive glasses as bone graft substitute. Välimäki VV, Aro HT. Scandinavian Journal of Surgery. 2006;95(2):95-102.
Osteoblast differentiation of bone marrow stromal cells cultured on silica gel and sol-gel-derived titania. Dieudonné SC, van den Dolder J, de Ruijter JE, Paldan H, Peltola T, van 't Hof MA, Happonen RP, Jansen JA. Biomaterials. 2002;23(14):3041-3051.
Histomorphometric and molecular biologic comparison of bioactive glass granules and autogenous bone grafts in augmentation of bone defect healing. Virolainen P, Heikkilä J, Yli-Urpo A, Vuorio E, Aro HT. J Biomed Mater Res. 1997;35A(1):9-17.

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






MRT

CT

© Helsinki University Central Hospital, Finland

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Literatur

Treatment of a recurrent aneurysmal bone cyst with bioactive glass in a child allows for good bone remodelling and growth. Lindfors NC. Bone. 2009;45:398-400.

Clinical experience on bioactive glass S53P4 in reconstructive surgery in the upper extremity showing bone remodelling, vascularization, cartilage repair and antibacterial properties of S53P4. Lindfors N. J Biotechnol Biomaterial. 2011;1(5) (An open access journal).

*nicht-osteinduktiv


© Department of Paediatric Orthopedic Surgery, Turku University Hospital, Finland

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Literatur

Treatment of a recurrent aneurysmal bone cyst with bioactive glass in a child allows for good bone remodelling and growth. Lindfors N. Bone, 2009;45(2):398-400.

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Röntgen 6 Monate postoperativ

Röntgen 6 Monate postoperativ


Courtesy of Dr. Ilkka Helenius, Turku University Hospital, Finland

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

Röntgen postoperativ

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© Department of Paediatric Orthopaedic Surgery, Turku University Hospital, Finland

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

Literatur

Bioactive glass S53P4 and autograft bone in treatment of depressed tibial plateau fractures. A prospective randomized 11-year follow-up. Perna K, Koski I, Mattila K, Gullichsen E, Heikkilä J, Aho AJ, Lindfors N. J Long-term Eff Med Impl. 2011;21(2):139-148.
Bioactive glass granules: a suitable bone substitute material in the operative treatment of depressed lateral tibial plateau fractures: a prospective, randomized 1 year follow-up study. Heikkilä JT, Kukkonen J, Aho AJ, Moisander S, Kyyrönen T, Mattila K. J Mater Sci Mater Med. 2011;22(4):1073-1080.
Posterolateral spondylodesis using bioactive glass S53P4 and autogenous bone in instrumented unstable lumbar spine burst fractures - A prospective 10-year follow-up study. Rantakokko J, Frantzén J, Heinänen J, Kajander S, Kotilainen E, Gullichsen E, Lindfors N. Scand J Surg. 2012; 101:66-71.

© Turku University Central Hospital, Finland

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






Röntgen 11 Jahre postoperativ

Röntgen 11 Jahre postoperativ

© Turku University Hospital, Finland

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Romanó et al. 2014

Literatur

A comparative study of the use of bioactive glass S53P4 and antibiotic-loaded calcium-based bone substitutes in the treatment of chronic osteomyelitis: a retrospective comparative study. Romanó CL, Logoluso N, Meani E, Romanó D, De Vecchi E, Vassena C, Drago L. Bone Joint J. 2014 Jun;96-B(6):845-850.
Bioactive glass BAG-S53P4 for the adjunctive treatment of chronic osteomyelitis of the long bones: an in vitro and prospective clinical study. Drago L, Romanó D, De Vecchi E, Vassena C, Logoluso N, Mattina R, Romanó CL. BMC Infectious Diseases 2013, 13:584 (An open access journal).
Through the looking glass; bioactive glass S53P4 (BonAlive®) in the treatment of chronic osteomyelitis. McAndrew J, Efrimescu C, Sheehan E, Niall D. Ir J Med Sci. 2013 Sep;182(3):509-511.
Bioactive glass S53P4 as bone graft substitute in treatment of osteomyelitis. Lindfors NC, Hyvönen P, Nyssönen M, Kirjavainen M, Kankare J, Gullichsen E, Salo J. Bone. 2010;47:212-218.

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Röntgen 1,5 Monate postoperativ

Röntgen 1 Jahr postoperativ

© Oulu University Hospital, Oulu Finland

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




1 Jahr postoperativ

1 Jahr postoperativ

© Infection Unit, Maastricht UMC+, the Netherlands

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
Röntgen CaP 2 Monate postoperativ

Röntgen BonAlive® granules 7 Monate postoperativ

© Helsingin University Central Hospital, Finland

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


Röntgen unmittelbar postoperativ

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
Röntgen unmittelbar postoperativ

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© Helsinki University Central Hospital

LITERATUR

Wirkmechanismus (Osteostimulation*)

Effects of bioactive glass S53P4 or beta-tricalcium phosphate and bone morphogenetic protein-2 and bone morphogenetic protein-7 on osteogenic differentiation of human adipose stem cells. Waselau M, Patrikoski M, Juntunen M, Kujala K, Kääriäinen M, Kuokkanen H, Sándor GK, Vapaavuori O, Suuronen R, Mannerström B, von Rechenberg B, Miettinen S. *J Tissue Eng.* 2012;3(1).

Osteoblast response to continuous phase macroporous scaffolds under static and dynamic culture conditions. Meretoja VV, Malin M, Seppälä JV, Närhi TO. *J Biomed Mater Res.* 2008;89A(2):317-325.

Molecular basis for action of bioactive glasses as bone graft substitute. Välimäki VV, Aro HT. *Scandinavian Journal of Surgery.* 2006;95(2):95-102.

Intact surface of bioactive glass S53P4 is resistant to osteoclastic activity. Wilson T, Parikka V, Holmbom J, Ylänen H, Penttinen R. *J Biomed Mater Res.* 2005;77A(1):67-74.

Granule size and composition of bioactive glasses affect osteoconduction in rabbit. Lindfors NC, Aho AJ. *J Mater Sci: Mater Med.* 2003;14(4):265-372.

Osteoblast differentiation of bone marrow stromal cells cultured on silica gel and sol-gel-derived titania. Dieudonné SC, van den Dolder J, de Ruijter JE, Paldan H, Peltola T, van 't Hof MA, Happonen RP, Jansen JA. *Biomaterials.* 2002;23(14):3041-3051.

Histomorphometric and molecular biologic comparison of bioactive glass granules and autogenous bone grafts in augmentation of bone defect healing. Virolainen P, Heikkilä J, Yli-Urpo A, Vuorio E, Aro HT. *J Biomed Mater Res.* 1997;35(1):9-17.

Bone formation in rabbit cancellous bone defects filled with bioactive glass granules. Heikkilä JT, Aho HJ, Yli-Urpo A, Happonen R, Aho AJ. *Acta Orthopaedica.* 1995;66(5):463-467.

Hemmung des Bakterienwachstums

***In vitro* antibiofilm activity of bioactive glass S53P4.** Drago L, Vassena C, Fenu S, De Vecchi E, Signori V, De Francesco R, Romanò CL. *Future Microbiol.* 2014;9(5):593-601.

Antibacterial effects and dissolution behavior of six bioactive glasses. Zhang D, Leppäranta O, Munukka E, Ylänen H, Viljanen MK, Eerola E, Hupa M, Hupa L. *J Biomed Mater Res.* 2010;93A(2):475-483.

Bactericidal effects of bioactive glasses on clinically important aerobic bacteria. Munukka E, Leppäranta O, Korkeamäki M, Vaahtio M, Peltola T, Zhang D, Hupa L, Ylänen H, Salonen JI, Viljanen MK, Eerola E. *J Mater Sci: Mater Med.* 2008;19(1):27-32.

Antibacterial effect of bioactive glasses on clinically important anaerobic bacteria *in vitro*. Leppäranta O, Vaahtio M, Peltola T, Zhang D, Hupa L, Ylänen H, Salonen JI, Viljanen MK, Eerola E. *J Mater Sci: Mater Med.* 2008;19(2):547-551.

***In situ* pH within particle beds of bioactive glasses.** Zhang D, Hupa M, Hupa L. *Acta Biomaterialia.* 2008;4(5):1498-1505.

Factors controlling antibacterial properties of bioactive glasses. Zhang D, Munukka E, Hupa L, Ylänen H, Viljanen MK, Hupa M. *Key Engineering Materials.* 2007;330-332:173-176.

Comparison of antibacterial effect on three bioactive glasses. Zhang D, Munukka E, Leppäranta O, Hupa L, Ylänen H, Salonen J, Eerola E, Viljanen MK, Hupa M. *Key Engineering Materials.* 2006;309-311:345-348.

Interactions between the bioactive glass S53P4 and the atrophic rhinitis-associated microorganism *Klebsiella ozaenae*. Stoor P, Söderling E, Grenman R. *J Biomed Mater Res.* 1999;48(6):869-874.

Antibacterial effects of a bioactive glass paste on oral micro-organisms. Stoor P, Söderling E, Salonen JI. *Acta Odontol Scand.* 1998;56(3):161-165.

Interactions between the frontal sinusitis-associated pathogen *Haemophilus Influenzae* and the bioactive glass S53P4. Stoor P, Söderling E, Andersson OH, Yli-Urpo A. *Bioceramics.* 1995;8:253-258.

*nicht-osteinduktiv

Chronische Osteomyelitis

A comparative study of the use of bioactive glass S53P4 and antibiotic-loaded calcium-based bone substitutes in the treatment of chronic osteomyelitis - a retrospective comparative study. Romanò CL, Logoluso N, Meani E, Romanò D, De Vecchi E, Vassena C, Drago L. Bone Joint J 2014;96-B:845-850.

Bioactive glass BAG-S53P4 for the adjunctive treatment of chronic osteomyelitis of the long bones: an *in vitro* and prospective clinical study. Drago L, Romanò D, De Vecchi E, Vassena C, Logoluso N, Mattina R, Romanò CL. BMC Infectious Diseases 2013;13:584. (An open access journal)

Through the looking glass; bioactive glass S53P4 (BonAlive®) in the treatment of chronic osteomyelitis. McAndrew J, Efrimescu C, Sheehan E, Niall D. Ir J Med Sci. 2013;182(3):509-511.

Clinical experience on bioactive glass S53P4 in reconstructive surgery in the upper extremity showing bone remodelling, vascularization, cartilage repair and antibacterial properties of S53P4. Lindfors NC. J Biotechnol Biomaterial. 2011;1(5). (An open access journal)

Bioactive glass S53P4 as bone graft substitute in treatment of osteomyelitis. Lindfors NC, Hyvönen P, Nyyssönen M, Kirjavainen M, Kankare J, Gullichsen E, Salo J. Bone. 2010;47:212-218.

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BonAlive Biomaterials Ltd
BonAlive Deutschland GmbH