TOUCH®

PLUS

new



Thanks to an increased amount of polyethylene, **the TOUCH® PLUS ensures enhanced durability under high daily demands**

CMC 1 PROSTHESIS

- **VITAMIN E**
  to limit polyethylene oxidation and maintain long-term mechanical performance
- **DUAL MOBILITY**
  to increase stability and range of motion
- **+35% POLYETHYLENE**
  compared with the TOUCH® STANDARD version



Fully committed to **Hand & Wrist** only

www.kerimedical.com

KeriMedical
medartis group

CMC 1 PROSTHESIS

TOUCH® PLUS



/ SPHERICAL CUPS

- 2 cup sizes (Ø 9 mm and 10 mm)
- 5 anti-rotation fins
- Spiked crown for stability
- Material: Stainless steel
- Dual coating: Titanium + HA



CTO09-PLUS



CTO10-PLUS

/ CONICAL CUPS

- 2 cup sizes (Ø 9 mm and 10 mm)
- Material: Stainless steel
- Dual coating: Titanium + HA



CTO109-PLUS



CTO110-PLUS

/ LINER AND STRAIGHT NECKS

- Liner material: Vitamin E Highly Cross-linked PE
- Neck material: Stainless Steel
- Straight neck
- 3 heights: 6, 8 and 10 mm



NTO06-PLUS



NTO08-PLUS



NTO010-PLUS

/ LINER AND OFFSET NECKS

- Liner material: Vitamin E Highly Cross-linked PE
- Neck material: Stainless Steel
- 15° offset neck
- 3 heights: 6, 8 and 10 mm



NTO156-PLUS



NTO158-PLUS




NTO1510-PLUS

/ STEMS

- 6 stems
- Anatomical trilobe section
- Material: Titanium
- Dual coating: Titanium + HA



STOXS STO0 STO1 STO2 STO3 STO4



/ INSTRUMENT SET

Complete set including orthopaedic standard instruments

TOUCH® PLUS rack



CE 2862

KERI MEDICAL SA
Chemin du Pré-Fleuri 5
1228 Plan-les-Ouates - Switzerland
Tel: +41 58 255 01 30
Fax: +41 58 255 01 40
customerservice@kerimedical.com

EU REP KERI MEDICAL FRANCE
34 rue Antoine Redier
74160 Archamps - France
Tel: +33(0)4 50 85 36 10
Fax: +33 (0)4 80 16 05 10
adv@kerimedical.com

Distributed by :

Sales brochure TOUCH® - SB-P_110-002.02_En - Issue date: 2026-06